Advantages

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Advanced Technology Delivers a 99.99% Sterility Assurance

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Intelligent Control The Brain of Equipment

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Energy & Eco-Friendly Enhancing Quality and Efficiency

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Ingenious Design Precision in Every Detail

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TIANS Equipment is committed to providing advanced energy-saving and environmentally friendly one-stop equipment solutions for customers in industries such as biological fermentation, biosynthetic, biopharmaceutical, food and health, biomaterials, and fine chemicals, etc.

Mastering multiple core technologies such as biological fermentation, sterilization, energy conservation and environmental protection, possessing dozens of energy-saving continuous sterilization systems and related patents for wastewater and waste gas treatment, and establishing a high standard and automated advanced equipment intelligent manufacturing center, our products are exported to multiple countries and regions at home and abroad, with a stable market share among the top.

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
Intelligent Instantaneous Sterilizing Equipment

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
设备 3D 演示
View 3D Video

TIANS Intelligent Instantaneous Sterilization Equipment works based on high-temperature sterilization technology, integrating high-efficiency heat exchangers, advanced holding tubes, and injectors, all operated through intelligent algorithms. With over 20 core patents, this equipment achieves the "Four Ultimate" standards — Ultimate Sterility, Ultimate Energy Efficiency, Ultimate Economic Benefit, and Ultimate Stability. As the national sales leader and recipient of the authoritative "Domestic Leading Level" certification, it is widely used for liquid material sterilization and inactivation of active substances in industries such as advanced materials, biopharmaceuticals, food and beverages.

Sterilization principle

According to the equivalent sterilization curve, the sterilization effect of conventional materials at 121°C, 30min is the same as that at 138°C, 36s. However, the destruction of nutrients in materials during disinfection is mainly determined by the sterilization time (except for materials that are extremely sensitive to heat).

According to this principle of high temperature and short time sterilization, after adding a certain safety factor, our company has designed this new type of continuous disinfection system by using the sterilization temperature of 138-140°C and the maintenance time of 50-60s.

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Classification & Composition

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■ 智能瞬时连消机分类方式 Classification method of intelligent instantaneous sterilizing equipment

根据加热方式不同分类 Classification according to different heating methods

- 直接加热
Direct heating
- 间接加热
Indirect heating

根据消毒方式不同分类 Classification according to different disinfection methods

- 蒸汽消毒 (空消)
Steam disinfection (steam disinfection)
- 过热水消毒 (水消)
Superheated water disinfection (water disinfection)

■ 智能瞬时连消机分为 4 个系列 Intelligent instantaneous sterilizing equipment is divided into 4 series

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Type I intelligent instantaneous sterilizing equipment - direct heating and steam disinfection

Type II intelligent instantaneous sterilizing equipment - direct heating and water disinfection

Type III intelligent instantaneous sterilizing equipment - indirect heating and steam disinfection

Type IV intelligent instantaneous sterilizing equipment - indirect heating and water disinfection

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Special runner size is selected, on-line cleaning is carried out, the structure is compact, no dead corners of materials exist, the resistance is small, no scaling, no coking matter and no material adhesion are realized.

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Make full use of latent heat of vaporization of steam, with compact structure, small resistance and no vibration.

Special flow channel size, high flow rate design and appropriate cleaning scheme can prevent the material from remaining due to the long-term operation of particles in the material.

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It is rolled from the whole plate, with no welding seam in the middle, and all automatic welding has no damage to the plate, and will not cause perforation or internal leakage of the plate.

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
Siemens S7-1200PLC control system is adopted and touch screen is equipped to facilitate on-site control. The intelligent system can multi-point sample the temperature, pressure and flow value of equipment with a sampling period of 100 ms. The system judges all temperature, pressure and flow signals and takes corresponding measures. For high-risk alarm, the system forcibly interrupts the process and switches to the safety mode. The system supports OPC/SQL standard database, can provide OPC, DP, MODBUS interface, can be used as SERVER end, or as CLIENT end. Necessary protocols and parameters are provided when necessary, which can be connected to the main DCS control system of the fermentation workshop for control.


Process

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
■ I 型智能瞬时连消机 - 直接加热空消智能瞬时连消机 Type I intelligent instantaneous sterilizing equipment


■ II 型智能瞬时连消机 - 直接加热水消智能瞬时连消机 Type II intelligent instantaneous sterilizing equipment

■ III 型智能瞬时连消机 - 间接加热空消智能瞬时连消机 Type III intelligent instantaneous sterilizing equipment

■ IV 型智能瞬时连消机 - 间接加热水消智能瞬时连消机 Type IV intelligent instantaneous sterilizing equipment

序号 No.	流量 (m³/h) Flow (m³/h)	进出口管径 (DN) Inlet and outlet pipe diameter (DN)	机组尺寸(长×宽×高 mm) Unit size (length × wide × height mm)	机组净重 (kg) Unit net weight (kg)	运行重量 (kg) Operating weight (kg)
1	3	32	3500×1700×2100	2500	2800
2	6	40	4000×2000×2100	3500	3900
3	10	50	4800×2100×2200	5000	5600
4	15	65	5100×2200×2200	5800	6800
5	30	80	5300×2500×3000	7800	9500
6	40	100	6200×2800×3000	12000	15000
7	50	100	6500×3000×3000	14000	17600
8	80	125	8500×3000×3000	24000	32500
9	100	150	9000×3000×3000	28000	40000

序号 No.	流量 (m³/h) Flow (m³/h)	进出口管径 (DN) Inlet and outlet pipe diameter (DN)	机组尺寸(长×宽×高 mm) Unit size (length × wide × height mm)	机组净重 (kg) Unit net weight (kg)	运行重量 (kg) Operating weight (kg)
1	3	32	5000×1700×2100	2700	3000
2	6	40	5500×2000×2100	3800	4200
3	10	50	6600×2100×2200	5500	6100
4	15	65	7000×2200×2200	6600	7600
5	30	80	7400×2500×3000	9400	11100
6	40	100	8500×2800×3000	14100	17100
7	50	100	8800×3000×3000	16600	20200
8	80	125	10800×3000×3000	28200	36700
9	100	150	11400×3000×3000	33200	45200